

Bleomycin: A Study of DNA Damage and the Cell Cycle

John Jeevarajan

Northwestern University




Dr. Honglu Wu

Dr. Ye Zhang

Radiation




SPACE LIFE SCIENCES
SUMMER INSTITUTE





Objectives of Internship


- Enhance proficiency in cell culture (human fibroblast cells)
- Become familiar with assay and analysis techniques
 - Immunofluorescence assays, fluorescence microscopy, DNA microarray analysis
- Apply learned techniques to further bleomycin ground study in preparation for parallel in-flight study
- Learn about space travel and science in space from SLSSI lectures



Background

- Space is a environment filled with radiation
 - Trapped protons, Galactic Cosmic Rays (GCRs), Solar Particle Events (SPEs)
- Assess astronaut's radiation risk
 - How cells respond to radiation
 - More important for long-term missions and missions outside of Near Earth Orbit (NEO)
- Question: How is DNA repair altered in microgravity compared to on Earth?
 - A controlled source of radiation is desired
 - Funded by NASA Fundamental Space Biology Program






The Space Radiation Environment



Space radiation : Energetic charged particles, high-LET (linear energy transfer)





Bleomycin

- Chemotherapy antibiotic
- Able to induce DNA double-strand breaks (DSBs) – radiomimetic
- Controlled source of DNA damage




γ -H2AX Assay and 53BP1 Assay

- Quantification of DSBs
 - Creation of foci as a result of DSBs from ionizing radiation
- H2AX – histone component variant
 - Phosphorylation of serine residue indicates DSB
 - Previous study tested high concentrations of bleomycin (0-80 $\mu\text{g/mL}$) and found little dose response (Venkatachalam et. al, 2008)
- 53BP1(p53-binding protein 1): protein involved in DNA repair and tumor suppression
 - Indirect quantification of DSBs




Cell Cycle Markers

- Determine if level of DNA damage correlates to stages of cell cycle
- Cyclin D1 and E: mid-to-late G_1 phase
- Cyclin A: high expression S phase, low expression in early G_2 phase
- Cyclin B1: G_2 phase. Expressed in nucleus in early M phase
- Phosphohistone H3: M phase
- Identify cell cycle stage that is best to study DNA damage and nuclear foci production



Bleomycin Treatment


- Growth of fibroblast cells in six 8-well tissue culture slides
- Once cells were 90% confluent, bleomycin treatment was performed
 - 0, 0.1, 1, and 10 $\mu\text{g/mL}$ for 3 hours
- Cells fixed with 4% paraformaldehyde
- Immunofluorescence staining:
 - 1: γ -H2AX and 53BP1
 - 2: γ -H2AX and Cyclin A (S phase)
 - 3: γ -H2AX and Cyclin B1 (G_2 phase)
 - 4: γ -H2AX and Phosphohistone H3 (M Phase)
 - 5: γ -H2AX and Cyclin D1 (G_1 phase)
 - 6: 53BP1 and Cyclin D1 (G_1 phase)



Results: γ -H2AX and 53BP1 Classification


Classification of cell types:

- Type I: Bright, pannuclear stain
- Type II: Uncountable bright foci or near pannuclear stain
- Type III: Countable foci
- Type IV: (only applies to 53BP1 analysis) γ -H2AX signal present, but observation of weak or no 53BP1 signal




Results: γ -H2AX and 53BP1 Foci Counts


- Distribution of foci counts: unimodal
- One type for all cell with countable foci



53BP1




γ -H2AX




53BP1 and γ -H2AX

γ -H2AX: Cell Distribution of Foci Counts




53BP1: Cell Distribution of Foci Counts



Results: γ -H2AX and 53BP1 Foci Counts

- Slight dose response for 53BP1, but not for γ -H2AX
- All significant changes except between 0.1 and 1 for γ -H2AX




Results: γ -H2AX and 53BP1 Foci Counts


□ As dose increases, Type I and II for γ -H2AX increase and Type I for 53BP1 increases




Type I



Type II




Type III




Type IV

γ -H2AX: % of Cells Per Type




53BP1: % of Cells Per Type




Results: Cell Cycle Study

- Cyclin D1 : G₁ phase
- Almost all cells are Type III
- Dose response observed for γ -H2AX foci


Cyclin D1: Dose vs. Cell Type Distribution




Cyclin D1: Dose vs. Average Number of Foci



Cyclin D1



γ -H2AX




Cyclin D1 and γ -H2AX



Results: Cell Cycle Study

- Cyclin A : S phase
- Most cells are Type III
- No dose response observed, due mainly to the high background in the control


Cyclin A: Dose vs. Cell Type Distribution




Cyclin A: Dose vs. Average Number of Foci



Cyclin A



γ-H2AX



Cyclin A and γ-H2AX


DNA Microarray

□ Analyzed microarray data from previous microarray samples

- GenePix4000B for data acquisition
- GenePix Pro 7 for data retrieval
- Microsoft Excel and EASE for analysis
- 0, 0.1, 1, and 10 ($\mu\text{g/mL}$)
- Up-regulated or down-regulated

□ DNA microarray system:

- Human GE 4x44K v2 Microarray
- ~44,000 transcripts
- Target 27,958 Entrez Gene RNAs



Significant Genes

Major categories affected (all concentrations):

- Apoptosis (cell death, programmed cell death, death)
- Regulation of Cell Proliferation
- DNA Repair or Response to DNA Damage/Endogenous Stimulus

EASE version 2.0 Academic use only

How to cite EASE | License terms and Disclaimer | Help | Exit

Note: before proceeding, you may want to ☒ update with the most recent online data

1. SELECT OUTPUT-> ☒ Display output in web browser
☐ Save the output to a file

2. INPUT GENES:
Either load
Affymetrix probesets
from a text file or paste below:
Clear all Paste
Count: 0

3. EXPLORE:
Either LINK genes to an online tool:
Link to: DAVID or link to NetAffx
or ANALYZE list:
Find Over-represented Gene Categories
☐ Enhance categorical information:
or ANNOTATE genes:
Annotate Genes
☐ Enhance annotation information:


Analysis options:

☐ Display links to PubMed evidence

Annotation options:

Refine Copy all
Make this a category of genes
Exit


Significant Genes




Significant Genes

- CDKN1A
 - Regulates G1 checkpoint by inhibiting activity of cyclin-CDK4 and cyclin-CDK2
 - Involved in p53 pathway in response to DNA damage
- MDM2
 - Part of autoregulatory negative feedback loop of the p53 pathway
 - Component in complex which links growth factor and DNA damage response pathways

CDKN1A Dose Response




MDM2 Dose Response




Significant Genes

- GADD45A
 - Increase in transcription usually follows environmental stresses and will lead to activation of the p38/JNK pathway
 - Activation of this gene due to DNA damage is mediated by p53 and other pathways
- DDB2
 - Required for DNA repair (especially UV-induced DNA repair)
 - Ubiquitination of histones H3 and H4 which aids cellular response to DNA damage

GADD45A Dose Response



DDB2 Dose Response



Conclusions

- 53BP1 and γ -H2AX are typically used markers for DNA DSBs caused by ionizing radiation, as indicated by colocalization of nuclear foci. However, similar to UV radiation, cells with bleomycin-induced DNA damage generally have more γ -H2AX foci than 53BP1 foci.
- 53BP1 and γ -H2AX have both shown dose dependent characteristics in identifying damage-induced foci. 53BP1 may be a better for DSBs because it showed a dose response more consistently.
- In terms of the cell cycle, cells in G_1 phase tend to have fewer damage-induced foci than cells in S phase.
- In general, an increase in number of genes involved in apoptosis, gene repair, and regulation of cellular processes is seen as the dose of bleomycin increases.
- Based on these data, we established that a concentration of 1 $\mu\text{g/mL}$ should be used for the flight study. In addition, cells in the G_1 phase should be analyzed for comparison to the ground study results.


Future Directions

- Conduct further gene studies based on DNA microarray data including pathway analysis and cluster analysis.
- Test bleomycin treatment for different lengths of time and allow specific time intervals for DNA repair to occur following treatment.
- Complete and analyze remaining cell cycle stains in order to conclude that ground and flight analysis should be completed based on G_1 phase cells.
- Become a practicing physician and aid the space industry in finding countermeasures to mitigate the effects of spaceflight on humans.



Acknowledgements

- Dr. Honglu Wu
- Dr. Ye Zhang
- Radiation Lab
- Judith Hayes
- Dr. Lauren Merkle
- Blythe Starkey
- Jackie Reeves
- Jan Connolly
- Dr. Ron McNeel
- Dr. Amanda Hackler



This work is partially funded by National Space Biomedical Research Institute via NASA Cooperative Agreement NCC 9-58

